

MCT PROTOCOLS



MCT SYSTEM

MCT Unit[®] | MCT Kit[®]

MCT Unit® and MCT Kit®




MCT Unit®

Device equipped to generate and apply electromagnetic and thermal energy across a broad range of temperatures and light wavelengths. Employing a photothermal biomodulation process, the MCT Unit® effectively improves existing regenerative medicine treatments. With multiple preset programs, the device is precisely tailored to the patient autologous material, resulting in the desired treatment outcome.

MCT Kit®

Patented cassette with a specific shape and chemical composition, exclusively developed for photothermal conditioning. Made of a medical grade polymer that guarantees optimal scattering, transmittance, and other optical properties, ensuring the emitted energy will reach the target effectively. Designed with a unique geometry that grants an excellent surface/volume energy exposure ratio, the MCT Kit® accommodates 10 mL of any autologous preparation.



MCT Unit® Specifications

The MCT Unit® is engineered for versatility, emitting light across a broad wavelength spectrum ranging from 467 nm to 850 nm. It is also equipped to thermally stimulate biological materials within a temperature range of 4°C to 45°C.



Photostimulation

- Max. emission energy 467 nm: 10 J/min.
- Max. emission energy 530 nm: 10 J/min.
- Max. emission energy 591 nm: 7.5 J/min.
- Max. emission energy 620 nm: 12 J/min.
- Max. emission energy 850 nm: 10 J/min.



Thermal stimulation


- Max. cooling temperature: 4°C
- Max. heating temperature: 45°C

This extensive scope of light and temperature options optimizes the flexibility and efficacy of the MCT system for different medical applications.

MCT Kit® Specifications

Open MCT Kit®

Closed MCT Kit®




The MCT Kit® has 2 syringe docking sites protected by 2 caps.

Docking sites are standard & Luer Lock® syringe compatible.

Insertion

Extraction




Insert material placing the MCT Kit® with caps facing upwards.

Extract material placing the MCT Kit® with caps facing downwards.

MCT Kit® Loading

The MCT Kit® is precisely crafted to fit in the tray of the MCT Unit® to enable proper closure.



MCT Unit[®] PRESETS

PRP | Cells | Exosomes

The MCT Unit® offers 3 distinct preset programs, each tailored to the specific characteristics of the autologous material and the desired final MCT product.



PRP



Cells



Exosomes

PRP



Goal

Stimulate platelet performance to achieve improved clinical results

Starting material

PRP

Program features

Stimulus 1

620 nm
1 J/cm²
10 min

Stimulus 2

4°C
15 min

Final material

MCT Plasma

The autologous material has to be **liquid** and contain **living structures**.

Cells



Goal

Promote cell metabolism for superior cellular-based therapy outcomes.

Starting material

Serum, biopsies, stromal vascular fraction, stem cells, and other cells.

Program features

Stimulus 1

620 nm
1 J/cm²
10 min

Stimulus 2

850 nm
0.2 J/cm²
10 min

Final material

MCT Cells

Exosomes



Goal

Induce exosome release, boosting the regenerative power of any product.

Starting material

Any product with living cells.

Program features

Stimulus 1

467 nm
2 J/cm²
10 min

Stimulus 2

37°C
10 min

Final material

MCT Exosomes



MCT PROTOCOL

MCT Plasma | MCT Cells | MCT Exosomes

MCT PROTOCOL

STEP 1



Obtain autologous material

- Procure autologous material.

STEP 2



Insert material into the MCT Kit®

- Position the MCT Kit® as in MCT Kit® Specifications page.
- Dock the syringe onto one docking site.
- Insert autologous material into the MCT Kit®.
- Close both caps of the MCT Kit®.

STEP 3




Set up the MCT Unit®

- Load the MCT Kit® in the MCT Unit® as in MCT Kit® Specifications page.
- Select the desired MCT preset:
[PRP](#) | [Cells](#) | [Exosomes](#)
- Press the Start button to start the session.
- During session, material undergoes a photothermal biomodulation process transforming it into **MCT material**.
- An acoustic alarm will sound upon completion of the session.
- Retrieve the MCT Kit® from the MCT Unit®.
- Open both caps of the MCT Kit® and dock the syringe as in MCT Kit® Specifications page.
- Extract the **MCT material** from the MCT Kit®.
Caution with spills from the other docking site.
- Administer the **MCT material** to the patient.

ADMINISTRATION TIMELINE

An interval of 3 weeks for treatment sessions and 4-6 months for maintenance sessions is advised, according to individual patient needs.






MCT CLINICAL PROTOCOL

PRP/Cells + Exosomes

PRP/Cells + Exosomes



Duration of complete treatment: **25–30 min**

1) Collect autologous material, aiming for a volume range of 18-20 mL. **2)** Split autologous material in two syringes (**syringe 1** and **syringe 2**), each containing 9-10 mL. **3)** Use **syringe 1** to fill the MCT Kit[®] with **half 1** of the material volume. **4)** Start PRP or Cells presets on the MCT Unit[®] and allow it to run upon completion of the session. Material undergoes a photothermal biomodulation process, transforming it into **MCT Plasma** (for PRP preset) or **MCT Cells** (for Cells preset). **5)** Retrieve **MCT material** from the MCT Kit[®] using **syringe 1**. **6)** Use **syringe 2** to refill the now empty MCT Kit[®] with the **half 2** of the material volume (9-10 mL). **7)** Start the Exosomes preset on the MCT Unit[®] and allow it to run upon completion of the session. **8)** While the MCT Unit[®] is running the Exosomes preset, administer to the patient the previously obtained **MCT Plasma** or **MCT Cells** material contained in **syringe 1**. **9)** Once the Exosomes preset is completed, autologous material has undergone photothermal biomodulation and the resulting material is now **MCT Exosomes**. Retrieve **MCT Exosomes** using **syringe 2**. **10)** Administer the **MCT Exosomes** material to the patient contained in **syringe 2**.



Inspiring technology
Outstanding results



M.C.T

Meta Cell Technology

Headquarters

Manel Farrés 101, 08173
Sant Cugat, Barcelona

Manufacturer

Manufactured by CLINIPRO, S.L.
Marketed by META CELL TECHNOLOGY, S.L.

Contact information

info@metacelltech.com
www.metacelltech.com



mctmetacelltech



metacelltech



+34 935 903 108

www.metacelltech.com

 **CLINIPRO SL**
Calle Manel Farrés, 101
08173 Sant Cugat del Vallès (Barcelona) - España